

BIOMÉRIEUX

# EMAG<sup>®</sup>

FLEXIBILITY. EFFICIENCY. CONFIDENCE.  
FOR ALL EXTRACTION NEEDS




PIONEERING DIAGNOSTICS


# HOW DO YOU MANAGE ALL YOUR EXTRACTION NEEDS EFFICIENTLY, WITH CONFIDENCE?

**EMAG® offers efficient and flexible extraction to manage both routine and challenging needs.**

**Count on bioMérieux's skilled team to give you the best support.**



## DIMENSIONS



## WEIGHT

EMAG®: 261 Kg  
Cabinet: 58 Kg

## UNMATCHED FLEXIBILITY

- 2 independent sections
- 1 set of reagents for all samples
- 1 unique disposable design

## INCREASED WORKFLOW EFFICIENCY

- 48 eluates in 98 minutes
- 5 predefined workflows
- 10 minutes hands-on-time

## CONFIDENCE & PEACE OF MIND

- Over 15 years of experience with > 250 publications
- > 5 safety parts against contamination
- 1 tailored service from installation to routine

## AUTOMATED NUCLEIC ACID EXTRACTION SYSTEM

### MARKING & COMPLIANCE

- CE
- IVD

### SAFETY


- UV lamp
- HEPA filter
- Liquid level detections

### CONNECTIVITY

- Bi-directionnal
- VILINK®

### EXPERTISE AND QUALITY

- BOOM technology
- Simplified mechanics



# EMAG<sup>®</sup>

## ORDERING INFORMATION



### INSTRUMENT

PRODUCT NAME	REFERENCE
EMAG <sup>®</sup>	Ref. 418591

### REAGENTS

PRODUCT NAME	REFERENCE
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Extraction Buffer 1	Ref. 280130
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Extraction Buffer 2	Ref. 280131
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Extraction Buffer 3	Ref. 280132
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Magnetic Silica	Ref. 280133
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Lysis Buffer	Ref. 280134

### CONSUMABLES

PRODUCT NAME	REFERENCE
NUCLISENS <sup>®</sup> easyMAG <sup>®</sup> Disposable	Ref. 280135
EMAG <sup>®</sup> Tips 1000 µL	Ref. 418922

Product availability varies by country.  
Please contact your local bioMérieux representative.